



BiomedicalInformatics


3184 Graves Hall 333 W. 10th Ave. Columbus, OH 43210 614-292-4778 (phone) 614-292-7659 (fax)

# Introduction to Microarray Data Analysis

- **Introduction to gene expression microarray**
  - A middle-man's approach
  - Applications of microarray
- **Microarray data processing/analysis workflow**
  - Data format and visualization
  - Data normalization
    - Two-color array
    - Affymetrix array
- **Software and databases**

# Review of Biology

mRNA, cDNA,  
exon, intron




# What is microarray?

- If we can assay every single molecule of DNA/RNA of interest directly, do we still need microarray?
- Currently direct single-molecule sequencing is still not mature, probes are used instead. Probe is a “middle-man”.


# How is microarray manufactured?

- Affymetrix GeneChip
  - silicon chip
  - oligonucleotide probes lithographically synthesized on the array
  - cRNA is used instead of cDNA




# How does microarray work?

RNA fragments with fluorescent tags from sample to be tested



# How does microarray work?




# **Two-major types of microarray**

- **Affymetrix-like arrays – single channel (background-green, foreground-red)**
- **cDNA arrays – two channel (red, green, yellow)**



# Affymetrix GeneChip

- silicon chip
- oligonucleotide probes lithographically synthesized on the array
- cRNA is used instead of cDNA





# Affymetrix GeneChip

- silicon chip
- oligonucleotide probes lithographically synthesized on the array




# Two-channel microarray

- Printed microarrays
- Long probe oligonucleotides (80-100) long are “printed” on the glass chip
- Comparative hybridization experiment



# Probe selection

- Protocol for extracting mRNA
- 3' bias – why? Think degradation.
- Multiple probes for one region







# How do we process microarray data (measurement)?

- cDNA array – ratio, log ratio

$$T_i = \frac{R_i}{G_i} \text{ OR } \log \text{ ratio} = \log_2 \frac{R_i}{G_i}$$



- Affymetrix array



$$\text{Difference}_{\text{probe pair}} = PM - MM$$

$$\text{Average Difference}_{\text{probe set}} = \sum_{i=1}^n \frac{(PM_i - MM_i)}{n}$$

# **Applications of microarrays**

- **Gene expression**
- **Exon expression**
- **SNP detection**
- **Copy number variance  
(arrayCGH)**
- **Tiling array (e.g., ChIP-chip)**

# Major vendors

- **Affymetrix**
- **Agilent**
- **Illumina**
- **Nimblegen**

- Introduction to gene expression microarray
  - A middle-man's approach
  - Applications of microarray
- **Microarray data processing/analysis workflow**
  - **Data format and visualization**
  - **Data normalization**
    - Two-color array
    - Affymetrix array
- Software and databases




# Typical workflow

- QC
- Normalization
- Visualization (boxplot, PCA, RI plot, etc).
- Comparative study (volcano plot)
- Clustering
- Network/pathway inference
- Motif finding

# Spatial Images of the Microarrays


- Data for the same brain voxel but for the untreated control mouse
- Background levels are much higher than those for the Parkinson's disease model mouse
- There appears to be something non random affecting the background of the green channel of this slide

F635 Median




$\times 10^4$

F532 Median




$\times 10^4$

B635 Median



10000  
8000  
6000  
4000  
2000

B532 Median



2500  
2000  
1500  
1000  
500


# Take a look ...




Poorly defined borders




Large holes




Saturated spot




Dust specs




Fiber or scratch?



Bubble



Edge effect





Background haze

(McShane, NCI)



# Example – Affymetrix Data Files

- Image file (.dat file)
- Probe results file (.cel file)
- Library file (.cdf, .gin files)
- Results file (.chp file)



# Example – Affymetrix Data Files

- Image file (.dat file)
- Probe results file (.cel file)





# Scatter plots of the Microarrays

- A measure of the **actual expression levels**, i.e., differences between the median foreground and the median background for the red channel and green channel:

"F635 Median - B635"


"F532 Median - B532"



# RI plots of the Microarrays

- RI (ratio-intensity) plot or MA plot



$$M = \log\left(\frac{R}{G}\right)$$



$$A = \frac{1}{2}(\log(R) + \log(G))$$

# Scatter plots of the Microarrays



## Bad Array Example




(McShane, NCI)



# Box plot



# Example:



# **Normalization – microarray data is highly noisy**

- **Intensity imbalance between RNA samples**
- **Affect all genes**
- **Not due to biology of samples, but due to technical reasons**
- **Reasons include difference in the settings of the photodetector voltage, imbalance in total amount of RNA in each sample, difference in uptaking of the dyes, etc.**
- **The objective is to adjust the gene expression values of all genes so that the ones that are not really differentially expressed have similar values across the**

# Two major issues to consider

- Which genes to use for normalization
- Which normalization algorithm to use

# Which genes to use for normalization

- Housekeeping genes
  - **Genes involved in essential activities of cell maintenance and survival**, but not in cell function and proliferation
  - These genes will be similarly expressed in all samples.
  - **Difficult to identify – need to be confirmed**
  - **Affymetrix GeneChip provides a set of house keeping genes** based on a large set of tests on different tissues and were found to have low variability in these samples (but still no guarantee).

# Which genes to use for normalization

- Spiked controls
  - Genes that are not usually found in the samples (both control and test sample). E.g., yeast gene in human tissue samples.

# Which genes to use for normalization

- Using all genes
  - Simplest approach – use all adequately expressed genes for normalization
  - **The assumption is that the majority of genes on the array are housekeeping genes and the proportion of over expressed genes is similar to that of the under expressed genes.**
  - If the genes on the chip are specially selected, then this method will not work.

# **Two-color array normalization**

- **Intra-slide normalization**
- **Inter-slide for cDNA arrays**





# Normalization

- **Linear (global) normalization**
  - Simplest but most consistent
  - Move the median to zero (slope 1 in scatter plot, this only changes the intersection)
  - No clear nonlinearity or slope in MA plot

$$X_i^{norm} = k * X_i$$



$$c = \log(k)$$

$$M_i^{norm} = \log(X_i^{norm}) = c + M_i$$




# Normalization

- **Intensity-based (Loess/Lowess) normalization**
  - Loess/Lowess fit
  - Overall magnitude of the spot intensity has an impact



# Normalization

- **Intensity-based normalization**
  - “Straighten” the Lo(w)ess fit line in MA plot to horizontal line and move it to zero





$$X_i^{norm} = k(A) * X_i$$

$$c(A) = \log(k(A))$$

$$M_i^{norm} = \log(X_i^{norm}) = c(A) + M_i$$

# Normalization



- **Intensity-based (Lowess) normalization**
  - Nonlinear
  - Gene-by-gene, could introduce bias
  - Use only when there is a compelling



(McShane, NCI)

# Normalization

- **Quantile normalization**
  - Nonlinear
  - Same intensity distribution





# Normalization

- Location-based normalization
  - Background subtracted ratios on the array may vary in a predictable manner.
  - Sample uniformly across the chip
  - Nonlinear
  - Gene-by-gene, could introduce bias
  - Use only when there is a compelling reason
- Other normalization method
  - Combination of location and intensity-based normalization

# Normalization

- Which normalization algorithm to use
  - Inter-slide normalization
  - Not just for Affymetrix arrays



# Normalization

- Linear (global) – the chips have equal median (or mean) intensity
- Intensity-based (Lowess) – the chips have equal medians (means) at all intensity values
- Quantile – the chips have identical intensity distribution
- Quantile is the “best” in term of normalizing the data to desired distribution, however it also changes the gene expression level individually
- Potential issues - overfitting



# **Affymetrix array normalization**

- **Inter-slide normalization only**
- **Probe-level normalization**
- **Affymetrix MicroArray Suite (MAS) 5.0**
- **Robust Multiarray Average (RMA)**
- **Quantile**
- **GC-RMA**

# Affymetrix array normalization

- **Inter-slide normalization only**
- **Probe-level normalization**
- **Affymetrix MicroArray Suite (MAS) 4.0**
  - Simple subtraction of MM from PM
  - Use only probes within 3 times of SD of PM-MM to exclude outliers
  - Not robust
- **MAS 5.0**
  - Use weight (Turkey Biweight Estimate) for each probe based on its intensity difference from the mean
  - Log transformed data for mean (geometric mean)
  - Robust

# Affymetrix array normalization

- **Robust Multiarray Average (RMA)**
  - Background correction on each chip.
    - Assuming strictly positive distribution. No negative numbers
    - Do NOT use MM information
  - Normalization (inter-chip).
    - Quantile
  - Probe level intensity calculation.
    - Linear model for signal, affinity, and noise.
  - Probe set summarization.
    - Combine probes for one probeset into a single number
    - Median polishing (chip to its median, gene to its median, iterate and converge)

# **Affymetrix array normalization**

- **GC-Robust Multiarray Average (GC-RMA)**
  - Correct back ground noise and non-specific binding
  - Affinity computed from position specific base effect
  - MM information is used (subtracted from PM after correction)

# **Affymetrix array normalization**

- **RMA/GCRMA pros and cons (comparing to MAS5.0)**
  - Less variance at low expression values
  - Less false positives
  - Consistent fold change estimates
  - More false negatives, especially for low-expression level probes
  - Quality control after normalization is difficult
  - Quantile normalization may overfit and hide real differences

- Introduction to gene expression microarray
  - A middle-man's approach
  - Applications of microarray
- Microarray data processing/analysis workflow
  - Data format and visualization
  - Data normalization
    - Two-color array
    - Affymetrix array
- **Software and databases**

# Microarray analysis software

- **Open source R**
- **Bioconductor**
- **BRBArray tools (NCI biometric research branch)**
- **Matlab Bioinformatics Toolbox**
- **Affymetrix Expression Console**
- **DChip**
- **GeneSpring**
- **Partek**
- **...**

# Microarray Databases

- **Gene Expression Omnibus (GEO) database – NCBI**
  - <http://www.ncbi.nlm.nih.gov/entrez/query.fcgi?DB=pubmed>
- **EMBL-EBI microarray database (ArrayExpress)**
  - <http://www.ebi.ac.uk/Databases/microarray.html>
- **Stanford Microarray Database (SMD)**
  - <http://genome-www5.stanford.edu/>
- **caARRAY sites**
- **The Cancer Genome Atlas (TCGA)**
- **Other specialized, regional and aggregated databases**
  - <http://psi081.ba.ars.usda.gov/SGMD/>
  - <http://www.oncomine.org/main/index.jsp>
  - [http://ihome.cuhk.edu.hk/~b400559/arraysoft\\_public.html](http://ihome.cuhk.edu.hk/~b400559/arraysoft_public.html)
  - ...



# Gene Expression Omnibus (GEO)

- <http://www.ncbi.nlm.nih.gov/projects/geo/query/browse.cgi>

## Total holdings

	Public	Unreleased	Total
Platforms	2727	319	3046
Samples	103186	24641	127827
Series	4351	980	5331

## Browse public holdings

- All contacts
- All platforms
  - in situ oligonucleotide (553)
  - spotted oligonucleotide (697)
  - spotted DNA/cDNA (1369)
  - antibody (5)
  - tissue (0)
  - MS (7)
  - SARST (1)
  - MPSS (7)
  - RT-PCR (6)
  - oligonucleotide beads (15)
  - mixed spotted oligonucleotide/cDNA (3)
  - spotted protein (1)
  - SAGE (38)
- All samples
  - RNA (94534)
  - genomic (6671)
  - protein (423)
  - SAGE (837)
  - mixed (403)
- All series

Oct.  
200  
6

## Total holdings

	Public	Unreleased	Total
Platforms	7925	517	8442
Samples	485908	87829	573737
Series	19157	3514	22671

## Browse public holdings

- All contacts
- All platforms
  - in situ oligonucleotide (2725)
  - spotted oligonucleotide (2021)
  - spotted DNA/cDNA (2445)
  - antibody (9)
  - tissue (0)
  - MS (16)
  - SARST (2)
  - MPSS (17)
  - RT-PCR (41)
  - oligonucleotide beads (128)
  - mixed spotted oligonucleotide/cDNA (12)
  - spotted protein (20)
  - SAGE (77)
- All samples
  - RNA (390777)
  - genomic (80401)
  - protein (2159)
  - SAGE (1707)
  - mixed (2259)
  - SRA (5454)
- All series

Oct.  
201  
0

# Gene Expression Omnibus (GEO)

- **GEO Profiles**

This database stores individual gene expression and molecular abundance profiles assembled from the [Gene Expression Omnibus \(GEO\)](#) repository. Search for specific profiles of interest based on gene annotation or pre-computed profile characteristics. GEO Profiles facilitates powerful searching and linking to additional information sources.

- **GEO DataSets**

This database stores curated gene expression and molecular abundance DataSets assembled from the [Gene Expression Omnibus \(GEO\) repository](#). Enter search terms to locate experiments of interest. DataSet records contain additional resources including cluster tools and differential expression queries.

(From GEO website)

# Gene Expression Omnibus (GEO)

- **GPL**

- A Platform record describes the list of elements on the array (e.g., cDNAs, oligonucleotide probesets, ORFs, antibodies) or the list of elements that may be detected and quantified in that experiment (e.g., SAGE tags, peptides). Each Platform record is assigned a unique and stable GEO accession number (GPLxxx). A Platform may reference many Samples that have been submitted by multiple submitters.

- **GSM**

- A Sample record describes the conditions under which an individual Sample was handled, the manipulations it underwent, and the abundance measurement of each element derived from it. Each Sample record is assigned a unique and stable GEO accession number (GSMxxx). A Sample entity must reference only one Platform and may be included in multiple Series.

# Gene Expression Omnibus (GEO)



- **GSE**

- A Series record defines a set of related Samples considered to be part of a group, how the Samples are related, and if and how they are ordered. A Series provides a focal point and description of the experiment as a whole. Series records may also contain tables describing extracted data, summary conclusions, or analyses. Each Series record is assigned a unique and stable GEO accession number (GSExxx).

- **GDS**

- GEO DataSets (GDS) are curated sets of GEO Sample data. A GDS record represents a collection of biologically and statistically comparable GEO Samples and forms the basis of GEO's suite of [data display and analysis tools](#). Samples within a GDS refer to the same Platform, that is, they share a common set of probe elements. Value measurements for each Sample within a GDS are assumed to be calculated in an equivalent manner, that is, considerations such as background processing and normalization are consistent across the dataset. Information reflecting experimental design is provided through GDS

# GEO Datasets



My NCBI  
[\[Sign In\]](#) [\[Register\]](#)

All DatabasesPubMedNucleotideProteinGenomeStructureOMIMPMCTournalsBooks

Search GEO DataSets for breast cancer Go Clear [Save Search](#)

Limits Preview/Index History Clipboard Details

Display Summary Show 20 Sort By  Send to

All: 781 DataSets: 89 Platforms: 27 Series: 665

Items 1 - 20 of 781 Page 1 of 40 Next

☐ **1: GDS3716 record: Breast cancer: histologically normal breast epithelium [ *Homo sapiens* ]** GEO Profiles, Links

Summary: Analysis of histological normal breast epithelia from both ER- and ER+ breast cancer patients and prophylactic mastectomy patients, and normal breast epithelia from reduction mammoplasty patients. Results provide insight into the mechanisms underlying breast cancer initiation and progression.  
Parent Platform: [GPL96](#)  
Reference Series: [GSE20437](#)


Type: Expression profiling by array, count

Subsets: 2 disease state, 4 specimen sets.

Supplementary Files: CEL [download...](#)

Samples: 42

- [GSM512539](#): reduction mammoplasty breast epithelium sample 1
- [GSM512540](#): reduction mammoplasty breast epithelium sample 2
- [GSM512541](#): reduction mammoplasty breast epithelium sample 3
- [GSM512542](#): reduction mammoplasty breast epithelium sample 4
- [GSM512543](#): reduction mammoplasty breast epithelium sample 5
- [GSM512544](#): reduction mammoplasty breast epithelium sample 6





☐ **2: GDS3638 record: Actin effect on breast cancer cell line: dose response and time course [ *Homo sapiens* ]** GEO Profiles, Links

▼ **Top Organisms [Tree]**

- Homo sapiens (669)**
- Mus musculus (113)
- Rattus norvegicus (18)
- Canis lupus familiaris (2)
- Macaca mulatta (1)
- Macaca fascicularis (1)
- Saccharomyces cerevisiae (1)
- Caulobacter vibrioides (1)


Recent activity


Turn Off Clear


-  [breast cancer](#) (781)
-  [gpl570](#) (1771) Geo Datasets

» See more...

# GEO Datasets



  
Gene Expression Omnibus



Search for

DataSet Record GDS3716:

Title:

Summary:

Organism:

Platform:

Citation:

Breast cancer: histologically normal breast epithelium

Analysis of histological normal breast epithelia from both ER- and ER+ breast cancer patients and prophylactic mastectomy patients, and normal breast epithelia from reduction mammoplasty patients. Results provide insight into the mechanisms underlying breast cancer initiation and progression.

*Homo sapiens*

GPL96: [HG-U133A] Affymetrix Human Genome U133A Array

Graham K, de las Morenas A, Tripathi A, King C et al. Gene expression in histologically normal epithelium from breast cancer patients and from cancer-free prophylactic mastectomy patients shares a similar profile. *Br J Cancer* 2010 Apr 13;102(8):1284-93. PMID: 20197764

Reference Series:

GSE20437

Sample count:

42


Value type:

count

Series published:

2010/06/01

Cluster Analysis



Download

DataSet SOFT file

Series family SOFT file

Series family MINiML file

Annotation SOFT file

Data Analysis Tools

Find genes ?

Compare 2 sets of samples

Cluster heatmaps

Experiment design and value distribution

Find gene name or symbol:

Find genes that are up/down for this condition(s):

☒ specimen

☒ disease state

NLM NIH GEO Help Disclaimer Section 508

# GEO Profiles

h GEO Profiles   [Save Search](#)

[Limits](#) [Preview/Index](#) [History](#) [Clipboard](#) [Details](#)

Display  Show  Subgroup effect

**All: 54675**

Items 1 - 20 of 54675

☐ 1: [GDS2250 record](#) | [GPL570](#) 204320\_at [Homo sapiens]  
Annotation: [COL11A1](#): collagen, type XI, alpha 1 CO11A1, COL16, STL2  
Reporter: [NM\\_001854](#)  
Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 2: [GDS2250 record](#) | [GPL570](#) 37892\_at [Homo sapiens]  
Annotation: [COL11A1](#): collagen, type XI, alpha 1 CO11A1, COL16, STL2  
Reporter: [NM\\_001854](#) [NM\\_080629](#) [NM\\_080630](#)  
Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 3: [GDS2250 record](#) | [GPL570](#) 204822\_at [Homo sapiens]  
Annotation: [TTK](#): TTK protein kinase ESK, FLJ38280, MPS1L1, PYT  
Reporter: [NM\\_003318](#)  
Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

# GEO Profiles

**Limits** **Preview/Index** **History** **Clipboard** **Details**

Display: Summary Show: 20 Subgroup effect Send to

All: 54675

Items 1 - 20 of 54675

Page 1 of 2734 Next

☐ 1: [GDS2250 record](#) | [GPL570](#) 204320\_at [Homo sapiens]  
Annotation: [COL11A1](#): collagen, type XI, alpha 1 CO11A1, COL6, STL2  
Reporter: [NM\\_001854](#)  
Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 2: [GDS2250 record](#) | [GPL570](#) 37892\_at [Homo sapiens]  
Annotation: [COL11A1](#): collagen, type XI, alpha 1 CO11A1, COL6, STL2  
Reporter: [NM\\_001854](#) [NM\\_080629](#) [NM\\_080630](#)  
Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 3: [GDS2250 record](#) | [GPL570](#) 204822\_at [Homo sapiens]  
Annotation: [TTK](#): TTK protein kinase ESK, FLJ38280, MPS1L1, PYT  
Reporter: [NM\\_003318](#)  
Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count


☐ 4: [GDS2250 record](#) | [GPL570](#) 205392\_s\_at [Homo sapiens]  
Annotation: [CCL14](#): chemokine (C-C motif) ligand 14 CC-1, CC-3, CKb1, HCC-1, HCC-3, MCIF, NCC-2, NCC2, SCYA14, SCYL2, SY14  
Reporter: [NM\\_004166](#)  
Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 5: [GDS2250 record](#) | [GPL570](#) 204719\_at [Homo sapiens]  
Annotation: [ABCA8](#): ATP-binding cassette, sub-family A (ABC1), member 8 KIAA0822


Internet



# GEO Profiles



- Left y-axis is (supposed to be) log two based (must check to verify) expression level.
- Right y-axis is the percentile of this expression level in the entire chip.
- All the chips are normalized.



# GEO Profiles

ases PubMed Nucleotide Protein Genome Structure PMC Journals Books

for GDS2250[ACCN] mmp14   [Save Search](#)

**Limits** Preview/Index History Clipboard Details

Display Summary  20

**All: 4**

Items 1 - 4 of 4

One page.

☐ 1: [GDS2250 record](#) | [GPL570](#) 160020\_at [Homo sapiens] 47 samples [Links](#)

Annotation: [MMP14](#): matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: [NM\\_004995](#)

Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 2: [GDS2250 record](#) | [GPL570](#) 202828\_s\_at [Homo sapiens] 47 samples [Links](#)

Annotation: [MMP14](#): matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: [NM\\_004995](#)

Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 3: [GDS2250 record](#) | [GPL570](#) 202827\_s\_at [Homo sapiens] 47 samples [Links](#)

Annotation: [MMP14](#): matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: [NM\\_004995](#)


Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

☐ 4: [GDS2250 record](#) | [GPL570](#) 217279\_x\_at [Homo sapiens] 47 samples [Links](#)

Annotation: [MMP14](#): matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: [NM\\_004995](#)

Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count



# GEO Profiles

ases PubMed Nucleotide Protein Genome Structure PMC Journals Books

for GDS2250[ACCN] mmp14 Go Clear Save Search

Limits Preview/Index History Clipboard Details

Display Summary Show 20 Mean Value Send to

All: 4

Items 1 - 4 of 4

One page.

47 samples Links

1: GDS2250 record | GPL570 160020\_at [Homo sapiens]

Annotation: MMP14: matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: NM\_004995

Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

2: GDS2250 record | GPL570 202828\_s\_at [Homo sapiens]

Annotation: MMP14: matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: NM\_004995

Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

3: GDS2250 record | GPL570 202827\_s\_at [Homo sapiens]

Annotation: MMP14: matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: NM\_004995


Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count

4: GDS2250 record | GPL570 217279\_x\_at [Homo sapiens]

Annotation: MMP14: matrix metalloproteinase 14 (membrane-inserted) MMP-X1, MT1-MMP, MTMMP1

Reporter: NM\_004995

Experiment: Basal-like breast cancer tumors, gene expression array-based transformed count




- Multiple probesets for different genes
- The number of probesets are different
- Probesets may have different versions
- May corresponding to polymorphism (splice variants)
- The results from different probesets may be inconsistent
- Various ways of combining the data

# **GEO Profiles**

- **Most new datasets are deposited as GSE series datasets instead of GDS datasets and cannot be visualized directly.**
- **Users need to download them for further processing.**
- **A simple way is to download the Data Matrix.**

# How do we use microarray?

- Profiling
- Comparative study
- Clustering
- Network inference



Supplementary Figure 1: Clustering of laboratory/platform combinations using log ratio values of common genes

# **Moving beyond microarray?**

- **Cut the middleman**
  - **Next generation sequencing**
  - **Single-molecule sequencing**
- **Where will microarray go?**
  - **Diagnosis**
  - **Specialized quick testing kit**